对折方式：M式对折（标题在上方）
纸张大小：95x145mm
打印方式：正反面打印

User Manual

22.04.22

【 Product Name】 Dental Implant Screwdriver

【 Model 】 C-SK1

【 Product Specifications 】

A100、B100、A119、B119、 A120、B120、A125、 B125、 A127、 B127、A140、 B140、 A155、B155、A170、B170.

【 Intended 】

1. The mandrel is used in dental implant surgery and its reatration, and during the preparation of dental surgery (enamel and dentine treatment , cutting crowns.etc).

2. The mandrel must only used in dental clinics or hospitals by professional medical personnel.

【 Contraindications 】

Those who have severe allergic reaction to stainless steel are prohibited.

【 Precautions 】

1. Keep the mandrel out of the reach of chidren .

2. Prohibit the use of hydrogen peroxide for cleaning or disinfection.

3. Please select the matching mandrel according to actual needs .

4. The corrosion resistance meets the requirements of 5.4 test evaluation class b in YY/T0149-2006.

5. The high-temperature detoxification performance requires that the surface and structure of the product have not changed after Pre-vacuum pressure steam sterilization .

【 Instruction 】

Put the mandrel into wrenth(torque)instument, and tighten or unscrew it with an appropriate force according to the product characteristics.

【 sterilization 】

1. Pre-vacuum pressure steam sterilization: Steriliza at 134 ℃, sterilize for 5 minutes, anddry for 30minutes ;

2. Down exhaust pressure steam sterilization method: Sterilizae at 132 ℃,steriliza for 20minutes, and dry for 40 minutes.

【 Maintenance 】

1. keep dry and avoid sunlight ;

2. After use, wash and fully dry and store indoors ;

3. Before use, check that the mandrel is not demaged or that no part of it is missing otherwise contact dealer or manufacturer.

【 Storage And Transportation Conditions 】

1. Storage temperature -10℃ - +55℃, relative humidity ≤90%

2. Keep it dry and straight up.

【 Guarantee 】

Mandrel is covered by 12 months warranty.

【 Symbols 】

	Warning / Caution		Consult instruction for use
	Products comply with WEEE directive		Keep dry
	Type B applied part		Thermal cleaning
	This way up	134°C	Sterilizable in at 134℃
	Manufacturer		

Foshan COXO Medical Instrument Co.,Ltd

BLDG 4, District A Guangdong New Light Source Industrial Base, South of Luocun Avenue Nanhai District Foshan 528226 Guangdong China

Tel:0757- 66692058 Web:www.coxotec.com

Version:Ver1.0 Date: 2022.04.14